

TOE/TEE kit (Adult)

REF 1250-01

Description

TOE/TEE kit non sterile.
A transducer cover (square tip) on an applicator, an extended 10 ml syringe filled with ultrasound transmission gel, a twist-lock and an adult bite guard.

Packaging information

Packaging	Peel pouch
Selling unit	02x24 kits
EAN13 (x24)	8717371750701
EAN13 (x48)	8717371750718

Specifications

Sterile	No
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	1100 mm



Probetection™ TOE/TEE kit.

Specially designed for Adult Trans Oesophageal/Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE kit (Adult)

REF 1250-02

Description

TOE/TEE kit sterile.
A transducer cover (square tip) on an applicator, an extended 10 ml syringe filled with ultrasound transmission gel, a twist-lock and an adult bite guard.

Packaging information

Packaging	Peel pouch
Selling unit	02x24 kits
EAN13 (x24)	8717371750725
EAN13 (x48)	8717371750732

Specifications

Sterile	Yes (gamma)
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	1100 mm



Probetection™ TOE/TEE kit.

Specially designed for Adult Trans Oesophageal /Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE cover (Adult)

REF 1250-03

Description

TOE/TEE cover Non sterile.
A single transducer cover
and a twist lock.
(Square tip)

Packaging information

Packaging	Sealed PE bag
Selling unit	05x50 covers
EAN13 (x50)	8717371750756
EAN13 (x250)	8717371750763

Specifications

Sterile	No
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	1100 mm



Probetection™ TOE/TEE cover

Specially designed for Adult Trans Oesophageal /Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE cover (Adult)

REF 1250-04

Description

TOE/TEE cover sterile.
A single transducer cover
and a twist lock.
(Square tip)

Packaging information

Packaging	Peel pack
Selling unit	05x24 covers
EAN13	NA

Specifications

Sterile	Yes (gamma)
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	1100 mm



Probetection™ TOE/TEE cover

Specially designed for Adult Trans Oesophageal /Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

Self foaming wash cloth REF 830

Description

Puracloth™, the rinse free, self foaming disposable wash cloth. Puracloth™ is a soft foam cloth for full and partial body cleansing and moisturizing. Puracloth™ is skin friendly, non-sensitizing and free of preservatives.

Packaging information

Packaging	Peel pack
PCS per pack	30 cloths
Selling unit	30 Packs
EAN13	871737174017

Specifications

Sterile	No
Material cloths	Polyurethane Foam
Latex free	Yes



Puracloth™ wash cloths.

Puracloth™, the rinse free, self foaming disposable wash cloth. Puracloth™ is a soft foam cloth for full and partial body cleansing and moisturizing. Puracloth™ is skin friendly, non-sensitizing and free of preservatives.

Instructions for use:

1: Open the packaging and remove a Puracloth™. 2: Moisten each individual Puracloth™ in warm or cold water. 3: Squeeze the Puracloth™ gently to activate the foam. 4: Wash only one area of the body to reduce the risk of cross contamination.

5: Discard the Puracloth™ in an appropriate way and repeat with a new cloth. We advise you to use 6 to 8 cloths for a full body wash. If feces is present, make sure to remove bulk prior to use. Do not use additional cleansers. Rinsing the skin is not required

Puracloth™ Features

- Neutral odour
- Time saving
- Efficient
- Cost effective
- A safe alternative to textile wash-mitts


Puracloth™

Purasoft™

Purasoft™ is a soft wash mitt ideal for a patient's daily wash. The tear proof, non-woven flannel, gives Purasoft™ superior softness, water absorption capacity and strength.

Purasoft™ is packed in a handy re-closable bag.

Because Purasoft™ is disposed of after each use, it reduces the risk of contamination during cleansing.

Purasoft™ features

- Time saving
- Efficient
- Cost effective
- A safe alternative to re-usable textile wash-mitts



Purasoft™ dry mitts

The safe alternative for textile mitts



Purawash™ impregnated mitts

Purawash™

Purawash™ is a soft, skin friendly impregnated disposable wash mitt for bathing and cleaning of immobile patients and patients with special needs. It is used when daily washing or cleaning of the patient is not possible due to immobility of the patient. Purawash™ is used in Acute Care, Long Term Care and Homecare environments.

Purawash™ features

- Time saving
- Cost effective
- A safe alternative to textile wash-mitts
- Impregnated with a skin friendly, PH neutral lotion



Puracloth™ wash cloths

Puracloth™

Puracloth™, the rinse free, self foaming disposable wash cloth. Puracloth™ is a soft foam cloth for full and partial body cleansing and moisturizing. Puracloth™ is skin friendly, non-sensitizing and free of preservatives.

Puracloth™ features

- Neutral odour
- Time saving
- Efficient
- Cost effective
- A safe alternative to textile wash-mitts

Trans-vaginal/rectal

REF 1\$00 series

Description

Trans-vaginal and trans-rectal transducers are widely used in gynaecological and urological procedures. This technique brings the transducer into closer proximity to the target organs providing optimal images for use in diagnosis.

Packaging information

Packaging Single packed
Selling unit \$) x%\$\$

Specifications

Sterile Bc
Non sterile Yes
Material Polyisoprene



Probetection™ trans-vaginal covers

The Probetection™ transducer covers offer an optimal fit for all trans-vaginal transducers currently in use. All transducer covers are individually packed, available non-sterile.

For these particular gynaecology procedures Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ trans-rectal covers

The use of trans-rectal transducers typically requires two different covers. One of these is a sanitary cover, the other is a water stand-off cover for near-field visualization.

The Probetection™ line of rectal transducer covers offers both types to facilitate the correct use of transducers in this application.

Trans-vaginal/rectal

REF 1100 series

Description

Trans-vaginal and trans-rectal transducers are widely used in gynaecological and urological procedures. This technique brings the transducer into closer proximity to the target organs providing optimal images for use in diagnosis.

Packaging information

Packaging Single packed
 Selling unit 24x144 (1100)
 24x100 (1102/6)
 Selling unit 05x24 (Sterile)

Specifications

Sterile Yes
 Non sterile Yes
 Material Latex



Probetection™ trans-vaginal covers

The Probetection™ transducer covers offer an optimal fit for all trans-vaginal transducers currently in use. All transducer covers are individually packed, available sterile and non-sterile.

For these particular gynaecology procedures Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ trans-rectal covers

The use of trans-rectal transducers typically requires two different covers. One of these is a sanitary cover, the other is a water stand-off cover for near-field visualization.

The Probetection™ line of rectal transducer covers offers both types to facilitate the correct use of transducers in this application.



Probetection™ latex non sterile ultrasound transducer covers.**REF 1100**

Cover rolled. Single packed. Ø Diameter Tip: 33 mm x 195 mm in length.

REF 1102

Cover rolled. Single packed. Ø Diameter Tip: 20 mm x 200 mm in length.

REF 1103

Cover rolled. Single packed. Ø Diameter Tip: 26 mm x 200 mm in length.

REF 1104

Cover rolled. Single packed. Ø Diameter Tip: 26 mm x 300 mm in length.

REF 1106

Cover rolled. Single packed. Ø Diameter Tip: 30 mm x 300 mm in length.

Latex Trans-vaginal/rectal covers

Probetection™ latex sterile ultrasound transducer covers**REF 1100-02**

Cover rolled. Single packed in a peel pack. Ø Diameter Tip: 33 mm x 195 mm in length.

REF 1102-02

Cover rolled. Single packed in a peel pack. Ø Diameter Tip: 20 mm x 200 mm in length.

REF 1103-02

Cover rolled. Single packed in a peel pack. Ø Diameter Tip: 26 mm x 200 mm in length.

REF 1104-02

Cover rolled. Single packed in a peel pack. Ø Diameter Tip: 26 mm x 300 mm in length.

REF 1106-2

Cover rolled. Single packed in a peel pack. Ø Diameter Tip: 30 mm x 300 mm in length.



Gyneacology/Urology

Intraoperative kit

REF 1237/9 series

Description

Specifically designed for draping phased arrays, linear arrays, and curved linear arrays during surgical procedures. These covers are packed with or without a gel syringe and two latex free elastic bands.

Packaging information

Packaging Single packed
Selling unit x24

Specifications

Sterile Yes/No
Material PE
Latex free Yes



Probetection™ intraoperative

Specifically designed for draping phased arrays, linear arrays, and curved linear arrays during surgical procedures. These drapes cover the probe head and cable, up to the system connector. These covers are single packed optionally including a gel syringe and two latex free elastic bands.

Use:

These covers are used during Biopsy, Abdominal Surgery, Cardiac Surgery and Vascular Surgery.

Probetection™ latex free Intraoperative kits

REF 1237-02

Intraoperative cover Sterile. Single packed in a peel pack including a short gel syringe and 2 elastic bands. Tip: 100 mm (flat) over 610 mm in length. Total length: 610 mm. Telescopically folded. (Box end)



REF 1239-02

Latex free Intraoperative kits

REF 1238-02

Intraoperative cover Sterile. Single packed in a peel pack including a short gel syringe and 2 elastic bands. Tip: 150 mm (flat) over 2440 mm in length. Total length: 2440 mm. Telescopically folded. (Box end)

REF 1239-02

Intraoperative cover Sterile. Single packed in a peel pack including a short gel syringe and 2 elastic bands. Tip: 150 mm (flat) over 1220 mm in length. Total length: 1220 mm. Telescopically folded. (Box end)



REF 1238-02



REF 1237, 38 & 39 serie

TOE/TEE kit (Adult)

REF 1251-01

Description

TOE/TEE kit non sterile.
A transducer cover (round tip) on an applicator, an extended 10 ml syringe filled with ultrasound transmission gel, a twist-lock and an adult bite guard.

Packaging information

Packaging	Peel pouch
Selling unit	02x24 kits
EAN13 (x24)	8717371750565
EAN13 (x48)	8717371750572

Specifications

Sterile	No
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	1100 mm



Probetection™ TOE/TEE kit.

Specially designed for Adult Trans Oesophageal/Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE kit (Adult)

REF 1251-02

Description

TOE/TEE kit sterile.
A transducer cover (round tip) on an applicator, an extended 10 ml syringe filled with ultrasound transmission gel, a twist-lock and an adult bite guard.

Packaging information

Packaging	Peel pouch
Selling unit	02x24 kits
EAN13 (x24)	8717371750589
EAN13 (x48)	8717371750596

Specifications

Sterile	Yes (gamma)
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	1100 mm



Probetection™ TOE/TEE kit.

Specially designed for Adult Trans Oesophageal /Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE cover (Adult)

REF 1251-03

Description

TOE/TEE cover Non sterile.
A single transducer cover
and a twist lock.
(Round tip)

Packaging information

Packaging	Sealed PE bag
Selling unit	05x50 covers
EAN13 (x50)	8717371750640
EAN13 (x250)	8717371750619

Specifications

Sterile	No
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	1100 mm



Probetection™ TOE/TEE cover

Specially designed for Adult Trans Oesophageal /Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE cover (Adult)

REF 1251-04

Description

TOE/TEE cover sterile.
A single transducer cover
and a twist lock.
(Round tip)

Packaging information

Packaging	Peel pack
Selling unit	05x24 covers
EAN13	NA

Specifications

Sterile	Yes (gamma)
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	1100 mm



Probetection™ TOE/TEE cover

Specially designed for Adult Trans Oesophageal /Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE cover (Adult)

REF 1251-05

Description

TOE/TEE cover Non sterile.
A single transducer cover on an applicator and a twist lock. (Round tip)

Packaging information

Packaging	Sealed PE bag
Selling unit	50 covers
EAN13 (x50)	8717371750640

Specifications

Sterile	No
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	1100 mm



Probetection™ TOE/TEE cover

Specially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE covers are available as a cover on an applicator and a twist lock.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE cover (Adult)

REF 1251-06

Description

TOE/TEE cover sterile.
A single transducer cover on an applicator and a twist lock. (Round tip)

Packaging information

Packaging	Peel pack
Selling unit	40 covers
EAN13 (x01)	8717371750725
EAN13 (x40)	8717371750831

Specifications

Sterile	Yes (gamma)
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	1100 mm



Probetecation™ TOE/TEE cover

Specially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE/TEE).

The Probetecation™ TOE/TEE covers are available as a cover on an applicator and a twist lock.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetecation™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE kit (Adult)

REF 1251-10

Description

Extended TOE/TEE kit sterile. A transducer cover (round tip) on an applicator, an extended 10 ml syringe filled with ultrasound transmission gel, a twist-lock and an adult bite guard.


Packaging information

Packaging	Peel pouch
Selling unit	02x24 kits
EAN13 (x24)	8717371751272
EAN13 (x48)	8717371751289

Specifications

Sterile	Yes (gamma)
Material	Polyurethane
Latex free	Yes
Tip Ø	15 mm
Tapered to Ø	32 mm
Total length	2000 mm

EXTENDED



Probetection™ TOE/TEE kit.

Specially designed for Adult Trans Oesophageal/Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE kit (Paediatric)

REF 1260-01

Description

Paediatric TOE/TEE kit Non sterile. A transducer cover (round tip) on an applicator, an extended 10 ml syringe filled with ultrasound transmission gel, a twist-lock and a paediatric bite guard.

Packaging information

Packaging	Peel pouch
Selling unit	02x24 kits
EAN13 (x24)	NA
EAN13 (x48)	NA

Specifications

Sterile	No
Material	Polyurethane
Latex free	Yes
Tip Ø	11 mm
Tapered to Ø	31 mm
Total length	1000 mm



Probetection™ TOE/TEE kit.

Specially designed for paediatric and small adult Trans Oesophageal /Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE kit (Paediatric)

REF 1260-02

Description

Paediatric TOE/TEE kit sterile. A transducer cover (round tip) on an applicator, an extended 10 ml syringe filled with ultrasound transmission gel, a twist-lock and a paediatric bite guard.

Packaging information

Packaging	Peel pouch
Selling unit	02x24 kits
EAN13 (x24)	NA
EAN13 (x48)	NA

Specifications

Sterile	Yes (gamma)
Material	Polyurethane
Latex free	Yes
Tip Ø	11 mm
Tapered to Ø	31 mm
Total length	1000 mm



Probetection™ TOE/TEE kit.

Specially designed for paediatric and small adult Trans Oesophageal /Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE kit (Neonatal)

REF 1280-01

Description

Neonatal TOE/TEE kit Non sterile.
A transducer cover (round tip) on an applicator, an extended 10 ml syringe filled with ultrasound transmission gel, a twist-lock and a paediatric bite guard.

Packaging information

Packaging	Peel pouch
Selling unit	02x24 kits
EAN13 (x24)	8717371750930
EAN13 (x48)	8717371750947

Specifications

Sterile	No
Material	Polyurethane
Latex free	Yes
Tip Ø	9 mm
Tapered to Ø	31 mm
Total length	980 mm



Probetection™ TOE/TEE kit.

Specially designed for neonatal and Paediatric Trans Oesophageal /Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

The 1280 range is suitable for all MICRO TOE/TEE transducer brands used. Amongst others the Philips S8-3t micro TEE Hitachi /Aloka UST-52110S and the Hitachi /Aloka UST-52120S

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.



Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard

TOE/TEE kit (Neonatal)

REF 1280-02

Description

Neonatal TOE/TEE kit sterile. A transducer cover (round tip) on an applicator, an extended 10 ml syringe filled with ultrasound transmission gel, a twist-lock and a paediatric bite guard.

Packaging information

Packaging	Peel pouch
Selling unit	02x24 kits
EAN13 (x24)	NA
EAN13 (x48)	NA

Specifications

Sterile	Yes (Gamma)
Material	Polyurethane
Latex free	Yes
Tip Ø	9 mm
Tapered to Ø	31 mm
Total length	980 mm



Probetection™ TOE/TEE kit.

Specially designed for neonatal and Paediatric Trans Oesophageal /Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

The 1280 range is suitable for all MICRO TOE/TEE transducer brands used. Amongst others the Philips S8-3t micro TEE Hitachi /Aloka UST-52110S and the Hitachi /Aloka UST-52120S

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Especially designed for Adult trans Oesophageal or Trans Esophageal Echocardiography (TOE / TEE)

The Probetection™ TOE/TEE is available as a complete ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE / TEE covers are made out of a latex free polyurethane material. All our TOE / TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

Puracloth™

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard